

Section 4 - Storage and handling

Samples are kept in the cold room and logged on arrival. The retention period is twelve months from the date of receipt unless the method states otherwise. Containers must be labelled with the sample number, the date, and the name of the person who took the sample.

Stock solutions are stable for twelve months when stored at TM -18 C. Check the seal before use and discard any container whose label is no longer legible. Record the check in the logbook.

Filters are supplied in two sizes. Use the larger one for turbid samples and the smaller one for everything else. Rinse the funnel with deionised water between samples to avoid carry-over.